

ارزیابی گازهای خونی به روش ساده، دقیق و قابل اطمینان



Available test configurations

	B-60	B	E	E-Cl	E-Ca	E-Glu	E-BUN	B-Lac
pH	X	X	X	X	X	X	X	X
pCO ₂	X	X	X	X	X	X	X	*
pO ₂	X	X	X	X	X	X	X	X
tHb		X	X	X	X	X	X	X
sO ₂		X	X	X	X	X	X	X
Na ⁺			X	X	X	X	X	
K ⁺			X	X	X	X	X	
Cl ⁻				X				
Ca ⁺⁺					X			
Glucose						X		
BUN (Urea)							X	
Lactate								X

*Pending FDA 510k clearance

OPTI CCA-TS2 Blood Gas Analyzer with single-use Cassettes

Portable Blood Gas Analyzer with Optical Fluorescence Technology

MEASURED PARAMETERS	RANGE	RESOLUTION	OPERATING PARAMETERS
pH	6.6 - 7.8 pH units	0.001 pH units	Minimum Sample Size: 125 µL, 60µL for B-60 cassette
pCO ₂	10 - 200 mmHg	0.1 mmHg	Sample Type: Whole Blood, Plasma, or Serum (depends on analyte tested)
pO ₂	10 - 700 mmHg	0.1 mmHg	Sample Application: Syringe, Capillary or ComfortSampler™
tHb	5 - 25 g/dL	0.1 g/dL	Sample Input: Automatic Aspiration
sO ₂	60% - 100%	0.1%	Analysis Time (Time to Results): < 120 Seconds
Na ⁺	100 - 180 mmol/L	0.1 mmol/L	Ambient Temperature: 10° C - 30° C (50°F - 86°F)
K ⁺	0.8 - 10 mmol/L	0.01mmol/L	Relative Humidity: 5% - 95% (Non-Condensing)
Ca ⁺⁺	0.2 - 3.0 mmol/L	0.01mmol/L	Available Languages: English, Spanish, French, Italian, German, Japanese, Chinese, Polish, Portuguese, Romanian, Russian, and Turkish
Cl ⁻	50 - 160 mmol/L	0.1mmol/L	
Glucose	1.7 - 22.2 mmol/L (30 - 400 mg/dL)	0.1mmol/L (0.1 mg/dL)	
BUN (Urea)	2.8 - 112.0 mmol/L (1 - 40 mmol/L)	0.1mmol/L (0.01 mg/dL)	
Lactate	0.3 - 17.5 mmol/L (2.7-157.7 mg/dL)	0.01 mmol/L (0.1 mg/dL)	
Barometric Pressure	300 - 800 mmHg	0.1 mmHg	
CALCULATED PARAMETERS	INPUT PARAMETERS		SECURITY
Actual Bicarbonate (HCO ₃ ⁻)	Operator ID	O ₂ Mode	Multiple levels of security with user-definable permissions
Standard Bicarbonate (st.HCO ₃ ⁻)	Patient ID	Vent Mode	QC Lockout: User Definable QC Lockout Functions
Base Excess (BE)	Accession Number	Pplat	Operator IDs: Storage for 300 operator IDs and passwords.
Base Excess Extra Cellular Fluid (BEecf)	Date of Birth	Mvol	
Base Excess Actual (BEact)	Sex	PIP	
Buffer Base (BB)	Temperature	Liter Flow	MISCELLANEOUS
Standard pH (st.pH)	Medical Record Number	Tvol	Size: 4.7 H x 14.2 W x 9.1 D in (12.0 x 36.2 x 23.0 cm)
Hydrogen Ion Conc. (cH ⁺)	Account Number	PS	Weight: 4.3 kg (9.6 lbs) with rechargeable Battery Pack
Total CO ₂ (tCO ₂)	Test ID	PEEP	Voltage/Frequency: 120 or 240 AC/50-60 Hz
Oxygen Saturation (SO ₂)	Patient Name	Rate	Power Consumption (Max.): 110 VA
Oxygen Content (O ₂ ct)	Age	CPAP	Interface: ASTM, ASCII, and POCT1 protocols (Ethernet and USB Types A/B)
Alveolar-arterial O ₂ Conc. (AaDO ₂)	Attending Physician	FIO ₂	Approvals: CE, UL, IVDD, FDA 510k cleared
PSO	Patient Location	MCHC	
St. Ionized Calcium (pH=7.4)(nCa ⁺⁺)	Sample Collection Time	RQ	
Hematocrit (Hct)	Sample Type	P50	
	Puncture Site	Bilevel Pressure	
	Allen's Test	I/E Ratio	
	tHb Type	Comment field	
	Bypass		



OPTI Medical Systems, Inc.
235 Hembree Park Drive
Roswell, Georgia 30076 USA
www.optimedical.com
+1-770-510-4444 or 800-490-6784

پارسه شرکت مهندسی پزشکی
تابان طب

تلفن: ۸۸۸۱۳۲۹۴-۵ فاکس: ۸۸۸۱۳۲۹۶